



Monthly Maintenance Checklist for Cardiac Science Automated External Defibrillators

Unit Serial #:	TFS Office:	POC:
----------------	-------------	------

Initial boxes as items are checked off. Keep completed form on file.

	Jan	Feb	Mar	Apr	May	June	July	Aug	Sept	Oct	Nov	Dec
<u>Step 1</u> Open Lid												
<u>Step 2</u> Status indicator should turn red												
<u>Step 3</u> Status indicator should turn back to green within 5 seconds												
<u>Step 4</u> Check expiration dates on pads												
<u>Step 5</u> Listen for voice prompts												
<u>Step 6</u> Close lid and confirm green status indicator												
Date												
Initials												

Daily Maintenance: Verify that the Status indicator is green.